

BT-550

Infant Warmer

- Far Infrared Heater
- 7inch LCD Display
- Apgar Timer
- Swivel Head
- LED Examination Lamp
- X-ray Tray

(OPTIONS)

- IR temperature measuring
- Lifting Stand
- Weighing Scale
- Masimo SpO2
- Wide Range of Accessories
- IV Plate
- Ringer Pole
- Basket



Technical Specification

Physical Characteristics

Size	• 1185 x 846 x 1877 mm
Mattress	• 495 x 810 x 27 mm
Weight	• 98Kg

Temperature Control

Control mode	• Pre-Warm mode/ Baby mode/ Manual mode
Skin temperature display range	• 20.0~45.0°C (±0.3°C)
Skin temperature control range	• 30.0~39.0°C (±0.3°C)
Skin temperature probe	• 2 available
Heater output setting range	• 0~100% (5% resolution)

Display

Display panel	• 7inch TFT color LCD
Alarm	• Visual and Audible Alarms
Apgar Timer	• 0~59min59sec, 1,5,10min Beep

Other features

LED Examination Lamp	• 10W *4, < 3000lm
Mattress platform tilting	• ±15°
Weighing Scale	• 0~10Kg, ±50g
Lifting Stand	• 592 ~ 792 mm
Basket	• 250 x 470 x 490 mm
SpO2 measurement range	• 1~100%
SpO2 accuracy	• 70~100 % : ±3 digits



7inch
TFT LCD
Display
with
Apgar Timer



LED
Examination
Lamp



IR Temperature
Measuring



Masimo
SpO2

Manufacturer: Bistos Co., Ltd.

7th Fl., A Bldg., 144-3, Sangdaewon-dong, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea.

Tel : +82 31 750 0340 (ext. 103)

Fax : +82 31 750 0344

E-mail : sales@bistos.co.kr

www.bistos.co.kr