







Portable design



SpO2



SpO2+NIBP



Trend Graphic

BT-720

- 4.3" color touch screen
- SpO2, Pulse, NIBP, Masimo SpO2
- Short/long trend graphic / trend table display
- Audio & visual alarm, adjustable alarm limit
- Automatic brightness adjustment
- Accurate SpO2 performance during motion and low perfusion
- Perfusion Index data / bar graph display
- Pitch tone variation for pulse rate
- Internal memory for data storage
- SD card for easy software upgradation
- Over 8 hours continuous working on rechargeable lithium-ion battery
- Specialized PC software for data review and analysis
- Option for Masimo SpO2, NIBP
- HL7 export to clinical information systems

Time	SpO2	PR	ET	FI	awRR
01-01-2016 00:00	---	---	---	---	---
01-01-2016 00:01	---	---	---	---	---
01-01-2016 00:02	---	---	---	---	---
01-01-2016 00:03	---	---	---	---	---
01-01-2016 00:04	---	---	---	---	---
01-01-2016 00:05	---	---	---	---	---
01-01-2016 00:06	---	---	---	---	---
01-01-2016 00:07	---	---	---	---	---
01-01-2016 00:08	---	---	---	---	---
01-01-2016 00:09	---	---	---	---	---
01-01-2016 00:10	---	---	---	---	---
01-01-2016 00:11	---	---	---	---	---
01-01-2016 00:12	---	---	---	---	---
01-01-2016 00:13	---	---	---	---	---
01-01-2016 00:14	---	---	---	---	---
01-01-2016 00:15	---	---	---	---	---
01-01-2016 00:16	---	---	---	---	---
01-01-2016 00:17	---	---	---	---	---
01-01-2016 00:18	---	---	---	---	---
01-01-2016 00:19	---	---	---	---	---
01-01-2016 00:20	---	---	---	---	---


Trend Tabular



PC Viewer Software

BISTOS Patient Monitor Technical Specification

Model		BT-720	
Category		Vital Sign Monitor	
Display		4.3" Touch LCD	
ECG	Lead Type		
	Gain Selection		
	Sweep Speed (mm/s)		
	Bandwidth : Diagnostic Mode		
	Monitoring Mode		
	Surgery Mode		
	Strong Mode		
	Heart Rate Range (bpm)		
Respiration	Method		
	Measurement Range		
	Sweep Speed (mm/s)		
SpO ₂	Measurement Range		0-100 %
	Accuracy (70-100%)	Adult/Pediatric	±2%
		Neonate	±3%
	Accuracy (0-69%)		unspecified
	Perfusion Index		0.05-20%
	Pulse Rate Range (bpm)		25-250
NIBP**	Method		Automatic Oscillometric
	Operation Mode		Manual / Auto / STAT
	Parameter		Systolic, Diastolic, Mean
	Systolic Range (mmHg)	Adult	40-270
		Pediatric	40-200
		Neonate	40-130
	Diastolic Range (mmHg)	Adult	10-210
		Pediatric	10-160
		Neonate	10-90
	Mean Range (mmHg)	Adult	20-230
		Pediatric	20-175
		Neonate	20-100
Temperature **	Range		0-50°C (41 to 122 °F)
	Parameter		T1, T2 and TD
IBP *	Channel		
	Range (mmHg)		
Printer *	Type		
	Print Speed (mm/s)		
	Paper size (mm)		
CO ₂ *	Sidestream		
	Mainstream		
Multi-gas/O ₂ *			
SpO ₂ -Masimo *		Masimo SpO ₂	
C.O. *	Method		
	Range		
Battery	Type (capacity)		Li-ion (4400mAh)
	Run Time		5 hour
	Charging Time		4 hour
PC Software Interface		RJ45, SD card slot	
Warranty		2 year	



BIO SIGNAL TOTAL SOLUTION

Bistos Co., Ltd. (Headquarter)

7th Fl., A Bldg., Woolim Lions Valley 5-cha, 302,
Galmachi-ro, Jungwon-gu, Seongnam-si,
Gyeonggi-do, Korea

Tel. +82 (31) 750 0340

Fax. +82 (31) 750 0344

bistos@bistos.co.kr

www.bistos.co.kr

Bistos America Inc. (US Branch)

22941, Triton Way, Suite #242, Laguna hills, CA., USA 92653

Tel. +1-949-614-8745

Fax. +1-949-614-8745

Site. www.bistosamerica.com

E-mail. info@bistosamericainc.com