

PORTABLE PHLEGM SUCTION UNIT

7A-23A



Over-flow
protection



Suction
liner



Oil-free
pump

- Input power: 180vA
- Max Vacuum: $\geq 0.06\text{MPa}$
- Max Air Flow: $\geq 20\text{L/min}$
- Noise: $\leq 60\text{dBA}$
- Operating time: 30 minutes on, 30 minutes off

SUCTION LINER

2L



PACKAGING PARAMETERS

Model	QTN (PCS/CTN)	G.W/N.W (KG)	Size (cm)
7A-23A	1	10.2/8.5	37.4×30.9×49.9

JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD

Manufacturing Address: No.1 Baisheng Road, Development Zone, Danyang,

Jiangsu Province, P.R.China

Tel: 86-0511-8690 0859 86-0511-8690 0860

E-mail: yuwell@yuyue.com.cn; sales@yuyue.com.cn

www.yuwell.com