

yuwell

You, We Love Life

Electric suction apparatus

7A-23B

Adopts completely plastic panel design, more fancy and fashionable, easy to carry, low noise, high negative pressure and large flux. It can be widely applied in the surgical operations which need to absorb phlegm, thick liquid etc.

Technical data:

Power Voltage: AC220V \pm 22V, 50Hz \pm 1Hz

Max negative pressure: \geq 90kPa

Operation noise: \leq 65 dB (A)

Input power: 180VA

Pumping frequency: ☐ 20 ☐ 20 L/min

Reservoir capacity: 2500mL/pc, 2pieces

G.W./N.W. : 22.5/19.5kgs

Packing : 1PC/CTN 45.5*40*86.5cm



Larger flow(20L)



Oil-free pump



Oil-free pump

Clean - no pollution of oil

Convenient - no need to add oil for maintenance

Safe - no positive pressure during using