



HK-400III

Syringe Pump

- Large & colorful LCD display
- History records
- RS232 interface
- Adjustable buzzer volume
- Anti-bolus function
- Various Visual and audible alarms
- 90° rotatable pole clamp convenient for horizontal bar & vertical IV pole
- Unique syringe intelligent recognition technology
- Dynamically display the pressure
- User self-define syringe brands
- Stackable function

HK-400III

Specification

Infusion mode	Rate mode, Time mode, Dose mode, intermittent infusion mode, TIVA Optional: Drug library
Applicable syringes	Various brands of sterile syringes:5ml,10ml,20ml,30ml,50ml/60ml
Syringe accuracy	±2%(after syringe correct calibration)
Infusion rate	5ml syringe:0.1–150ml/h 10ml syringe:0.1–300ml/h 20ml syringe:0.1–600ml/h 30ml syringe:0.1–900ml/h 50ml syringe:0.1–1500ml/h
Volume to be infused	0~9999.9(ml)
Volume Infused	0~9999.9(ml)
KVO rate	0.1~2ml/h adjustable;default value:0.1ml/h
Bolus rate	5ml syringe:100–150ml/h 10ml syringe:100–300ml/h 20ml syringe:100–600ml/h 30ml syringe:100–900ml/h 50ml syringe:100–1500ml/h
Alarms	Infusion completion,Empty,Occlusion,Near completion,Low battery, Syringe disengaged,No AC power supply,Battery exhausted,Malfunction etc
Power supply	DC:12±1.2V AC:100~240V,50/60Hz
Power consumption	25VA
Nurse call alarm	Optional
Classification	Class I,type CF
Dimensions	230*137*139 (Length x Width x Height)
Weight	1.7KG
Waterproof Level	1P×3
Occlusion pressure	(40~160)kPa
Battery	Li_polymer 7.4V 1900Ah, running about 6 hours at 5ml/h
Operating conditions	Temperature:5℃~40℃ Relative humidity:10~95% Atmosphere pressure:86~106kPa
Transport conditions	Temperature:–20℃~60℃ Relative humidity:10~95% Atmosphere pressure:50~106kPa
Storage conditions	Temperature:–20℃~45℃ Relative humidity:10~95% Atmosphere pressure:50~106kPa



SHENZHEN HAWK MEDICAL INSTRUMENT CO., LTD.

Address:2B,Building No.2,Aerospace Science & Industry factory areas,
Yousong community,Longhua street,Bao'an District,Shenzhen,China

E-mail : szhk@hawkmedical.cn
TEL : 0086-0755-83151901
FAX : 0086-0755-83151906
Web : www.hawkmedical.cn

