



The "LIFETIME" Surgical Aspirator has been designed to use in operation theaters during surgical operation to take off the accumulated liquids in the operated resign of the patient and for application in Gynecology and Dermatology (liposuction).

The electronically controlled foot paddle and selection of variable capacity accumulation jars gives wide range working facilities to the operator.

SPECIFICATION

MOTOR TYPE	OILLESS PISTON PUMP
VOLTAGE	230 VAC \pm %10 , 50 Hz
MAX. VACUUM (ADJUSTABLE)	- 550 mmHg \pm %10
MAX. FREE AIR FLOW RATE	12 LPM \pm %10
NOISE LEVEL	< 50 dB
JARS CAPACITY	1 lt / 2 lt OPTIONAL
CLASSIFICATION	ISO 10079-1

ACCESSORIES

