



MX-5000

Vital Signs Monitor

FEATURES

- Compact and flexible appearance
- Simple and convenient operation, which is applicable for community medical treatment and family ect.
- 4.3 inch color TFT display and highlight LED display Adjustable audible and visual alarms
- Built-in rechargeable lithium battery for uninterrupted monitoring

NIBP

- Method: Oscillometry
- Measure mode: Manual, Auto, STAT
- Unit: mmHg, kPa
- Measure and alarm range:

Adult Mode

SYS 40 ~ 280 mmHg

DIA 10 ~ 225 mmHg

MEAN 20 ~ 240 mmHg

- Pediatric Mode

SYS 40 ~ 220 mmHg

DIA 10 ~ 160 mmHg

MEAN 20 ~ 170 mmHg

- Neonatal Mode

SYS 40 ~ 135 mmHg

DIA 10 ~ 100 mmHg

MEAN 20 ~ 110 mmHg

- Resolution: 1mmHg
- Accuracy: $\pm 5\text{mmHg}$

SPO2

Item

Specification

- Measuring Range: 0 ~ 100 %
- Resolution: 1%
- Accuracy: 70% ~ 100% ($\pm 2\%$)
- Pulse Rate: 20-250 BPM
- Resolution: 1 BPM
- Accuracy: $\pm 3\text{BPM}$

TEMP

- Measure and Alarm Range: 0 ~ 50 °C
- Resolution: 0.1°C
- Accuracy: $\pm 0.1^\circ\text{C}$

STANDARD PARAMETERS

- NIBP, SPO2, PR, TEMP

Packing

- Single package size: 35 *25 *26cm
- Single gross weight: 3 KG
- 8 unit per carton, package size: 72.5 *53 *56.5 cm, total gross weight 24 KG

