



# MX-GY

The Multifunctional  
Electric Surgical & Delivery Bed  
Precision, Performance  
and Practicality



## FEATURES

- Adjustable table surface (back, hip, and leg plates)
- Smooth-lift column and stable base
- Efficient transmission system with intuitive handheld controller
- Practical accessories: arm rests, leg supports, and mounts for anesthesia infusion stands

## SPECIFICATIONS

- Table Length × Width: 1780 (±30) mm × 600 (±30) mm
- Table Height (Min.): 755 (±20) mm
- Table Tilt: Forward tilt  $\geq 20^\circ$  ; Backward tilt  $\geq 17^\circ$
- Head/Back Section (Integrated): Upward fold  $\geq 30^\circ$
- Leg Section: Extendable
- Arm Support: (550 ±10) mm × (160 ±10) mm
- Leg Support: (275 ±10) mm × (190 ±10) mm
- Power Supply: 220V, 50Hz
- Maximum Load Capacity: 270 kg

## INDEPENDENT POWER SYSTEMS

ach function—height adjustment, tilt, and backrest positioning—is driven by its own motor, ensuring safe, flexible, and stable operation. Sleek, easy-to-clean covers and waste basins enhance hygiene and aesthetics.

## PRODUCT HIGHLIGHTS

- Precision Engineering – Equipped with CE approved motors and electric actuators for exceptional rigidity, smooth movement, and accurate positioning. Quiet, durable, and built for long-term reliability.
- High-Grade Construction – The bed frame is crafted from premium nickel-chromium stainless steel, combining elegance with effortless cleaning and sterilization.
- Comfort & Durability – A leather-upholstered sponge tabletop and a one-piece molded mattress provide lasting comfort, style, and resilience.
- Clinically Optimized – Removable auxiliary panels allow seamless adaptation to procedures, supporting doctors and midwives in dynamic clinical environments.

